



Product Lines

Laboratory Chemical Analyzers



Full Automated Biochemistry Analyzer

****Other Analyzers are also Available**

Application:

- Blood Lipid
- Kidney Function
- Liver Function
- Total Blood Protein
- Cardiac Enzyme
- Electrolyte
- Brain-Spinal Fluid

Features:

- Rapid response testing
- Fully automated
- Online dilutions
- 220 tests/hour
- Control of the entire automated process
- 3 pipette tips for dispensing
- Unattached reagent and sample probe
- Liquid level detection
- Collision protection
- Micro sample size
- Stability, precision and accuracy
- Automatic washing system
- Automatic storage and backup

Swine Influenza A Virus (IVA H1N1) Diagnostic Kit

Features:

- 24 tests per kit (25uL reaction)
- Rapid detection of clinically relevant subtypes
- Pre-mixed PCR reagents
- Applicable to most commercially available Real-Time instruments
- Kits have same instrument set-up parameters



****Other PCR Reagent Kits are also Available**

High Performance Liquid Chromatograph (HPLC)

Features:

- Wide wavelength range UV-VIS spectrometer
- Isocratic pump
- Dual piston design
- Floating pump seal
- Self priming feature
- Intuitive software control
- Easy maintenance to all parts

Application:

- Drugs
- Biological and natural compounds
- Food analysis
- Environmental analysis



Smart Solutions for your Laboratory needs

KGC RESOURCES SDN BHD (223165-D)
 2-4-2, Jalan Setia Prima E, U13/E
 Seksyen U13, Setia Alam,
 40170 Shah Alam

Tel: 603-3341-2880
 Fax: 603-3343-9880
 Email: sales@kgcscientific.com
www.kgcscientific.com